

**PAJUNK®**


## *SonoSystem*

*The complete system for  
ultrasound guided nerve blocks*




Nerve blocks

## SonoLong Sono



Set consisting of:

- SonoLong NanoLine cannula with Luer Lock connector and electrical connecting cable
- Adaptable injection tube
- SonoLong Sono catheter 20 G with integrated stainless steel helical coil, steel stylet and central opening
- Clamping adapter (yellow)
- Bacterial filter 0.2 µm
- FixoLong



Product	Size	SonoLong Sono catheter	Item No.	PU
<b>SonoLong NanoLine cannula</b>				
Facet tip	19 G x 50 mm	20 G x 50 cm	531187-31A	10
	19 G x 75 mm	20 G x 50 cm	561187-31A	10
	19 G x 100 mm	20 G x 50 cm	521187-31A	10
	19 G x 150 mm	20 G x 50 cm	511187-31A	10
SPROTTE® SPECIAL tip	19 G x 60 mm	20 G x 50 cm	531187-31B	10
	19 G x 120 mm	20 G x 50 cm	521187-31B	10
Tuohy tip	18 G x 50 mm	20 G x 50 cm	531187-31C	10
	18 G x 75 mm	20 G x 50 cm	561187-31C	10
	18 G x 100 mm	20 G x 50 cm	521187-31C	10
	18 G x 150 mm	20 G x 50 cm	511187-31C	10

# Complete SonoSystem


## The systems at a glance

### StimuLong Sono II



Set consisting of:

- SonoLong NanoLine cannula with Luer Lock connector and electrical connecting cable
- Adaptable injection tube
- StimuLong Sono catheter 20 G x 50 cm with central opening, electrically conductive stylet and integrated metal helical coil
- StimuLong clamping adapter (yellow) with integrated stimulation connection
- Connecting cable
- Bacterial filter 0.2µm
- FixoLong



Product	Size	StimuLong catheter	Item No.	PU
<b>SonoLong NanoLine cannula</b>				
Facet tip				
	19 G x 50 mm	20 G x 50 cm	531187-32A	10
	19 G x 100 mm	20 G x 50 cm	521187-32A	10
Tuohy tip				
	18 G x 50 mm	20 G x 50 cm	531187-32C	10
	18 G x 100 mm	20 G x 50 cm	521187-32C	10